

# ONE TRAY<sup>®</sup>

Sealed Sterilization Container



ONETRAY.COM



# Why ONE TRAY®?

Simply put, we are more efficient and have a lifetime warranty.

**ONE TRAY®** is an American made, lifetime warranted sealed sterilization container offering efficient instrument reprocessing in a fraction of the time it takes sterile wrap or other rigid containers. Mitigate delays caused by sterile wrap tears or wet packs and experience time and cost savings across multiple departments within your facility. Don't think about how to incorporate **ONE TRAY®** into your daily schedule, think about how you can incorporate your daily schedule into **ONE TRAY®**.

**Efficient, Effective, Economical.**



# Tray Design

**ONE TRAY**®'s build quality is designed to minimize dents, nicks and other damage. Our rigid container consists of a base, lid, filter covers, and deck plate constructed of anodized aluminum and stainless steel materials.



Compression Force Cam Locking design provides optimal lid to base interface

Deck plate elevates medical devices off the bottom of the tray


Filter Covers use a patented **SterObex™** barrier providing sealing capabilities for sterility maintenance



*Up close view of **SterObex™** Barrier Technology*

Unique lip design helps prevent nicks to the top edge of the base

Base design provides stackability for easy storage after sterilization



- Patented design uses thermodynamics to dispense steam throughout the container
- Staggered vent placement allows steam to reach lateral edges quickly
- Single intake vent with dual exhaust vents creates effective vapor exchange performance




# Processing Kits

The **ONE TRAY®** processing kits include single use disposable ☼ filters that permits the penetration of steam during the sterilization cycle and serves as a bacterial and fluid barrier at the conclusion of the sterilization cycle.

Kits include:

- 3 Filters - disposable hydrophobic SMS Filters
- 1 Chemical Integrator
- 2 Locks - with indicator dots
- 2 Load Cards - for easy identification

*The validated performance of **ONE TRAY®** is in effect **ONLY** with the exclusive use of **ONE TRAY®** components, filters, and tamper evident locks. Medical Devices should be prepared and sterilized according to the device manufacturer's written instructions for sterilization exposure parameters.*



**ONE TRAY®**  
Sealed Sterilization Container

# Tray Models and Descriptions

Container Catalog #	M2104	M2106	M2108	M2404	M2406	M2408
Description	21"l x 12"w x 4"h with 21"deck plate	21"l x 12"w x 6"h with 21"deck plate	21"l x 12"w x 8"h with 21"deck plate	24"l x 12"w x 4"h with 24"deck plate	24"l x 12"w x 6"h with 24"deck plate	24"l x 12"w x 8"h with 24"deck plate

## Steam Sterilization Parameters

Parameter	Dynamic Air Removal/ Pre-Vacuum Steam Sterilization Cycle
Exposure Temperature	270°F (132°C)
Exposure Time	4 minutes
Dry Time	Not Applicable
Cycle Cool Time	Varies according to load
Stacking	Not permitted in sterilization cycle



Please reference the **ONE TRAY®** IFU for usage instructions.



## Questions About Our Technology?

**ONE TRAY®** has a team of industry and product experts that can answer any questions.

Join our **ONE TRAY® Questions** Facebook group | **937.619.0138** | [info@onetray.com](mailto:info@onetray.com)

*Our mission is to provide the highest quality sterilization products supported by years of reliable service to our customers. We are committed to educating our industry on how new technology affects current practices and guidelines to help create an environment that benefits not only the user but supports our customers in their effort to provide One Standard of Care.*

# ONE TRAY®

**Sealed Sterilization Container**



Rev01  
condoc\_336



Innovative Sterilization Technologies  
7625 Paragon Road, Suite A, Dayton, OH 45459  
937.619.0138 | [www.onetray.com](http://www.onetray.com)



Please reference the **ONE TRAY®** IFU for usage instructions.



Made in the USA