

ONE TRAY® MODELS AND SIZES

Original 510k Models		Current Product Codes
ONE 20MH	equivalent to	M2108
ONE 20MHL	equivalent to	M2106
ONE 20MLH	equivalent to	M2106
ONE 20ML	equivalent to	M2104
ONE 23MH	equivalent to	M2408
ONE 23MHL	equivalent to	M2406
ONE 23MLH	equivalent to	M2406
ONE 23ML	equivalent to	M2404

	Original 510k	Current Product Code	Original 510k	Current Product Code	Original 510k	Current Product Code	Original 510k	Current Product Code
Model	ONE 20ML	M2104	ONE 20MLH	M2106	ONE 20MHL	M2106	ONE 20MH	M2108
Height "	4.00	5.50	5.50	7.50	6.00	7.50	7.50	9.50
Width "	11.50	11.62	11.50	11.62	11.50	11.62	11.50	11.62
Length "	20.50	20.88	20.50	20.88	20.50	20.88	20.50	20.88

	Original 510k	Current Product Code	Original 510k	Current Product Code	Original 510k	Current Product Code	Original 510k	Current Product Code
Model	ONE 23ML	M2404		M2406	ONE 23MHL	M2406	ONE 23MH	M2408
Height "	4.00	5.50	5.50	7.50	6.00	7.50	7.50	9.50
Width "	11.50	11.62	11.50	11.62	11.50	11.62	11.50	11.62
Length "	23.50	23.88	23.50	23.88	23.50	23.88	23.50	23.88

#REF!