



ONE TRAY® Containers Physical and Dimensional Specifications

Model	Revised standardized product description based on the following OUTSIDE drawing dimensions	Base	Drawing Dimension Length	Drawing Dimension Height	Lid	Drawing Dimension Length	Drawing Dimension Height	Total Height	
M2104	(L) 21.253 x (H) 5.378	CASE-2125-02	20.877	3.391	CASE-2115-01	21.253	1.996	5.387	n/a
M2106	(L) 21.253 x (H) 7.387	CASE-2145-02	20.877	5.391	CASE-2115-01	21.253	1.996	7.387	n/a
M2108	(L) 21.253 x (H) 9.387	CASE-2145-02	20.877	5.391	CASE-2135-01	21.253	3.996	9.387	n/a
M2404	(L) 24.253 x (H) 5.387	CASE-2425-02	23.877	3.391	CASE-2415-01	24.253	1.996	5.387	n/a
M2406	(L) 24.253 x (H) 7.387	CASE-2445-02	23.877	5.391	CASE-2415-01	24.253	1.996	7.387	n/a
M2408	(L) 24.253 x (H) 9.387	CASE-2445-02	23.877	5.391	CASE-2435-01	24.253	3.996	9.387	n/a

Model	Revised standardized product description based on the following INSIDE drawing dimensions plus deck plate	Base	Drawing Dimension Length	Drawing Dimension Height	Lid	Drawing Dimension Length	Drawing Dimension Height	Total Height without Deck Plate	Total Height with Deck Plate
M2104	(L) 20.142 x (H) 4.573	CASE-2125-02	20.877	3.328	CASE-2115-01	20.142	1.933	5.261	4.573
M2106	(L) 20.142 x (H) 6.573	CASE-2145-02	20.877	5.328	CASE-2115-01	20.142	1.933	7.261	6.573
M2108	(L) 20.142 x (H) 8.573	CASE-2145-02	20.877	5.328	CASE-2135-01	20.142	3.933	9.261	8.573
M2404	(L) 23.142 x (H) 4.573	CASE-2425-02	23.877	3.328	CASE-2415-01	23.142	1.933	5.261	4.573
M2406	(L) 23.142 x (H) 6.573	CASE-2445-02	23.877	5.328	CASE-2415-01	23.142	1.933	7.261	6.573
M2408	(L) 23.142 x (H) 8.573	CASE-2445-02	23.877	5.328	CASE-2435-01	23.142	3.933	9.261	8.573

Filter Area Specs:

3 filter areas

Filter area per - 7.625" x 8.375"

Total square inches per - 63.85938

Total container filter area - 191.5781

